

Palm Bladder Scanner

Noninvasive, Accurate, Reliable, Easy to Use



SASYACHOOK

AUSTRALIA

Ph: 1300 411 973 Fax: 03 9877 1643
Email: sales@sasyachook.com.au
www.sasyachook.com.au

Specification

- ◆ Transducer freq: 2.6MHz
- ◆ Scan method: 3D Mechanical sector scanning
- ◆ Scan depth: $\geq 140\text{mm}$
- ◆ Sway angle: $120^{\circ}\pm 2^{\circ}$
- ◆ Rotating angle: $180^{\circ}\pm 2^{\circ}$
- ◆ Display screen: 8.0" Color LCD TFT screen
- ◆ Power supply: Battery
- ◆ Patient Case storage: 100pcs
- ◆ Patient ID: 12 digits
- ◆ Max. Measurement: 999ml
- ◆ Scan mode: male; female; PRV



Application

PBS is widely used in hospital
Urology, Obstetrics and Gynecology,
Anorectic
Rehabilitation, Emergency
Department, General Surgery
Pediatrics, ICU

PBS's latent application area:
Appertain area, Uric incontinence
guide, SCI/MS
Mentality healing institution, Rest
home
Family tend:
SCI/MS, Community Medical Center,
Family and Community Doctor.





Weight: 2.1kg (N/W)



Palm Bladder Scanner

Mianyang Meike Electronic Equipment Co.,Ltd

Add: No.238, Chuangye Road, kechuangyuan Zone, Mianyang City, Sichuan, China



Noninvasive, Accurate, Reliable, Easy to Use

SASYACHOOK

AUSTRALIA

Ph: 1300 411 973 Fax: 03 9877 1643
Email: sales@sasyachook.com.au
www.sasyachook.com.au