

IndeeLift®

The Human Floor Lift

For Professional Healthcare

*Acute Hospitals
Skilled Nursing
Assisted Living*

The People Picker Upper!
HFL-400-P/550-P



"IndeeLift is the tool you need to safely get your residents off the ground after a fall."





The People Picker Upper!

HFL-400-P/550-P

Purpose

The HFL-400-P/550-P models are portable human floor lifts designed for use in professional healthcare environments to safely lift patients and residents who have fallen and are unable to get up from the floor. These sturdy, stable lifts utilize medically certified components with a working loads of 181kg and 250kg respectively.

Operation

After a fall, the healthcare professional wheels the HFL to the fallen person and assists them onto the seat. With the simple push of a button on the wired remote, the fallen person is safely and easily raised from the floor.

Summary of Features

- Lifts a person weighing up to 181 / 250kg to a height of 53.4cm above the floor.
- No complex cradles or harness to configure.
- Easy maneuverability for flexible positioning in tight spaces including ADA-compliant bathrooms.
- Large ramp-mount seat for easy slide-on access.
- Swinging/removable Rise Assist Handles for easy transfer to wheelchair.
- Powder coated for easy cleaning and decontamination.
- Powered by a battery pack with built-in charger.
- Audible indicator when charging is required.
- Charges from any standard AC power outlet.



Product Specifications

Physical Specifications	HFL-400	HFL-550
Safe Working Load.....	400lbs.....181kg	550lbs250kg
Maximum Overall Length.....	25".....635mm	29"..... 735mm
Maximum Overall Height.....	62".....1575mm	64".... 1626mm
Minimum Overall Height.....	41".....1041mm	45".... 1143mm
Maximum Seat Height.....	21".....534mm	21".....534mm
Turning Radius.....	30".....762mm	32".....813mm
External Width	23".....584mm	26.25".....667mm
Overall Height of Legs.....	2.5".....64mm	2.5".....64mm
Rear Casters.....	3".....76mm	5".....127mm
Power Pack.....	6.6 lbs..... 3kg	6.6lbs.....3kg
Total Weight.....	96lbs.....44kg	108 lbs.....49kg

All measurements are within a +5/-5% tolerance.

Electric Shock Protection

Charger – Class II

Lift – Internal power source

Degree of Shock Protection

Charger – Type B

Lift – Type B

Operating environment: >+5°C <+40°C

Electrical Specifications

Batteries – 2x12 volt Rechargeable sealed lead acid type

Battery Capacity – 3.2 Amp Hours

Charger Rated Input – 100/240VAC 50/60Hz

Charger Rated Output – 27.4/29 VDC@0.8A

T: 1300 411 973

E:sales@sasyachook.com.au

